

TOTEM

WWW.ROTAID.COM



 **ROTAID™**
AED CABINETS

www.rotaid.com

OFFICIAL ACCESSORIES

 **ROTAID™**
AED CABINETS

ROTAID TOTEM

PRODUCT SHEET

The Rotaid AED Cabinets totem allows for a stand alone in- and outdoor installation of your defibrillator. The solid construction and high quality appearance make sure your defibrillator stands out and can be securely placed in practically all environments and situations.

The modular set-up ensures an easy installation and the totem can easily be demounted when used for training and trade show purposes.

On the front, back and top of the totem, information and branding panels can be mounted. These panels are custom made, connecting to the customer's request and location, making sure the totem aligns perfectly with its surroundings.

The cabinet is safely mounted on the totem and comes with a locking system. The solid aluminium construction allows for installation in the most demanding environments and makes sure the functioning and appearance is guaranteed for long-lasting use.

+

APPLIANCES

- Training and fares
- Stand alone installations
- Locations with no wall mounting possibilities
- Outdoor installations
- Renting and eventing

+

DESIGN

- Custom made front, back and top panel
- Full colour UV protected printing
- Easy click-on mounting
- Allows for fixed floor mounting
- In combination with Rotaid Plus, Solid Plus and Solid Plus Heat cabinets

- For indoor and outdoor usage
- Solid and safe construction
- Locked fixation of the cabinet
- Custom made designs and panels

DIMENSIONS	
Footplate	680 mm round - 20 mm high Material: Steel powdercoated
Uprights / Stairs	2 legs of 640 mm Two steps of 480 mm Material: anodized aluminum
Panels	290 x 1380 mm Material: dibond 3 mm thick (aluminum / PVC)



ROTAID TOTEM

TECHNICAL DETAILS

